


ns >>

AUTO 10mm/mV



10mm/mV

aVR



aVL



aVF



aVF



2.5mm/mV

V1



V2



V3



V3




5mm/mV



80mmx20m-30m


25mm/s AC50Hz+EMG25Hz



80mmx20m-30m



25mm/s



2024-05-28 15:13

ID :
Nom :
Sexe :
Taille : cm
SYS : mmHg
DIA : mmHg

Âge :
Poids : kg

FC [bpm] : 79
Intervalle PR [ms] : 199
Durée P [ms] : 120
Durée QRS [ms] : 81
Durée T [ms] : 167
QT/QTc [ms] : 338/388
Axe P/QRS/T [deg] : 66.4/81.9/24.0
R(V5)/S(V1) [mV] : 1.88/2.57
R(V5)+S(V1) [mV] : 4.45

<< Conclusion

Ryth. sinusoï
Axe électriqu

Médecin