

AUTO 10mm/mV



10mm/mV



5mm/mV




II 10mm/mV



80mmx20m-30m

25mm/s AC50Hz+EMG25Hz




80mmx20m-30m

10mm/mV



25mm/s

10mm/mV



10mm/mV



2024-05-18 17:10

ID :
Nom :
Sexe :
Taille: cm
SYS : mmHg
DIA : mmHg

Âge :
Poids: kg

FC [bpm] : 70
Intervalle PR [ms] : 145
Durée P [ms] : 120
Durée QRS [ms] : 73
Durée T [ms] : 173
QT/QTc [ms] : 329/355
Axe P/QRS/T [deg] : 73.0/-36.2/59.9
R(V5)/S(V1) [mV] : 0.63/1.14
R(V5)+S(V1) [mV] : 1.77

<< Conclusion

Ryth. sinusoï
Dévia. ax. gc

Médecin