


AUTO 10mm/mV



10mm/mV




10mm/mV



II 5mm/mV



25mm/s AC50Hz+EMG25Hz



2024-05-18 16:50

<< Conclusions >

ID :
 Nom :
 Sexe :
 Taille : cm
 SYS : mmHg
 DIA : mmHg

Âge :
 Poids : kg

Ryth. sinusoïdal
 Dévia. ax. gche

FC [bpm] : 91
 Intervalle PR [ms] : 191
 Durée P [ms] : 120
 Durée QRS [ms] : 70
 Durée T [ms] : 174
 QT/QTc [ms] : 341/419
 Axe P/QRS/T [deg] : 54.3/10.8/16.5
 R(V5)/S(V1) [mV] : 2.38/1.39
 R(V5)+S(V1) [mV] : 3.77

Médecin