

AUTO 10mm/mV



10mm/mV



10mm/mV



II, 10mm/mV



25mm/s AC50Hz+EMG25Hz

80mmx20m-30m

5mm/mV

V4

V5

V6



10mm/mV

I

II

III

aVR

25mm/s




10mm/mV

aVL

aVF

V1

V2



10mm/mV

V3

V4

V5

V6



2024-05-18 18:14

ID : _____

Nom : _____

Sexe : _____

Taille : _____ cm

SYS : _____ mmHg

DIA : _____ mmHg

Âge : _____

Poids : _____ kg

FC [bpm] : 76

Intervalle PR [ms] : 187

Durée P [ms] : 120

Durée QRS [ms] : 74

Durée T [ms] : 187

QT/QTc [ms] : 346/390

Axe P/QRS/T [deg] : 59.2/34.4/38.7

R(V5)/S(V1) [mV] : 1.49/0.68

R(V5)+S(V1) [mV] : 2.17

<< Conclusion

Ryth. sinuso

Axe électriq

Médecin