

AUTO 10mm/mV




10mm/mV



5mm/mV




II 10mm/mV



25mm/s AC50Hz+EMG25Hz

80mmx20m-30m

5mm/mV





10mm/mV



10mm/mV



5mm/mV



25mm/s

2024-05-18 16:14

<< Conclusions >>

ID : _____

Nom : _____

Sexe : _____

Taille : _____ cm

SYS : _____ mmHg

DIA : _____ mmHg

Âge : _____

Poids : _____ kg

FC [bpm] : 83

Intervalle PR [ms] : 172

Durée P [ms] : 120

Durée QRS [ms] : 64

Durée T [ms] : 177

QT/QTc [ms] : 326/382

Axe P/QRS/T [deg] : 74.6/56.0/42.6

R(V5)/S(V1) [mV] : 1.67/0.85

R(V5)+S(V1) [mV] : 2.52

Ryth. sinusoïdal n
Axe électrique cou

80mmx20m-Médecin