

Parameter	Pred.	Value1	%Pred. 1	Value2	%Pred. 2	Value3	%Pred. 3	Parameter
FVC(L)	4.22	3.03	71.71	2.79	66.05	2.90	68.64	FVC(L)
FEV0.5(L)	/	2.09	/	2.23	/	1.45	/	FEV0.5(L)
FEV0.5/FI	/	69.22	/	80.06	/	50.15	/	FEV0.5/FIVC(%)
FEV1(L)	3.62	2.79	77.09	2.74	75.88	2.50	69.07	FEV1(L)
FEV1/FIVC	70.92	92.13	129.92	98.45	138.83	86.24	121.61	FEV1/FIVC(%)
FEV3(L)	4.22	2.94	69.67	2.75	65.13	2.64	62.68	FEV3(L)
FEV3/FIVC	85.70	94.83	110.65	99.84	116.50	94.45	110.20	FEV3/FIVC(%)
FEV6(L/S)	/	97.16	/	98.61	/	91.31	/	PIF(L/S)
FEV0.2-1.	6.00	1.55	25.83	1.17	19.50	1.49	24.83	FIF0.2-1.2(L/S)
FEV50(L/S)	/	2.94	/	2.75	/	2.64	/	FIF50(L/S)
FEV25-75(L)	/	97.11	/	98.59	/	91.31	/	FIF25-75(L/S)
PEF(s)	8.28	6.11	73.79	7.04	84.93	3.46	41.82	FIT(s)
FEF0.2-1.2(L/S)	/	5.56	/	6.67	/	2.63	/	EVOL(ml)
FEF25(L/S)	7.02	6.11	87.12	7.04	100.25	2.65	37.69	MVV(L/M)
FEF50(L/S)	4.66	4.26	91.46	4.87	104.58	3.45	74.14	ELA(Year)
FEF75(L/S)	2.19	2.01	91.60	2.75	125.46	2.47	112.90	
FEF25-75(L/S)	4.35	3.79	87.15	4.69	107.84	3.11	71.47	
FEF75-85(L/S)	/	1.73	/	2.40	/	2.17	/	
MET(s)	/	0.33	/	0.22	/	0.71	/	

Clinical reference only

Chc 341

Pred.	Value1	%Pred. 1	Value2	%Pred. 2	Value3	%Pred. 3
4.22	0.18	4.22	0.85	20.22	2.68	63.60
/	0.08	/	0.25	/	1.41	/
/	44.04	/	28.94	/	52.45	/
3.62	0.14	3.94	0.46	12.65	2.31	63.86
70.92	79.88	112.64	53.62	75.62	86.05	121.34
/	0.18	/	0.84	/	2.34	/
/	99.56	/	98.97	/	87.38	/
8.28	0.18	2.17	0.61	7.38	3.68	44.38
/	/	/	/	/	2.00	/
/	0.13	/	0.43	/	3.53	/
/	0.14	/	0.37	/	3.33	/
/	1.51	/	3.09	/	1.81	/
/	87.38	/	39.16	/	251.60	/
126.52	97.56	77.11	96.03	75.90	87.42	69.09
21.00	66.00	314.29	68.00	323.81	75.00	357.14

